

# Evaluating Wheat and Corn Flours for Tortilla Production

Dr. M. Hikmet Boyacioglu  
Cereal Scientist  
Applications Development Specialist  
KPM Analytics  
hboyacioglu@kpmanalytics.com



TIA Tech Conference  
October 13-14, 2021  
Long Beach, CA, USA



# Introduction

---

## TORTILLA UPDATE

Value, flexibility, and health and wellness innovation spur category growth!

The value and flexibility of tortillas have proved to be a winning formula for the grain-based foods category during the past year, as the coronavirus pandemic-led consumers to search out more inexpensive and convenient ways to feed their families.

As a result, sales have surged.

The shift by tortilla manufacturers toward value-added products began several years ago, but the momentum remains apparent in recent sales figures.

Schroeder, 2021



# Introduction

---

## TORTILLAS

A tortilla is a flat, round, unfermented bread produced from wheat flour or lime-cooked maize!

Wheat flour tortillas are more prevalent than maize tortillas throughout the rest of the world out of Central America.

Processing and characteristics of **wheat** flour and **maize** tortillas differ considerably.

Rooney and Serna-Saldivar, 2016

# Introduction

## WHEAT VS CORN/MAIZE



	Wheat	Corn
Taxonomic family	Triticum	Zea
	Both are belonging to the grass family "Graminea (Poacea)" Grown for their edible seeds	
Origin	Mesopotamia	Mesoamerica
Products	Flour types, semolina, baked goods, pasta products, bulgur, ...	Corn starch, syrup, grits, meal, oil, flour, masa, Tortillas, snacks, ...
Production (20/21)	773 million tons	1,127 million tons
Constituents	High protein quantity and quality, Gluten!	High carbohydrate, oil



# Wheat Flour Tortillas

---

## TYPES OF COMMERCIALY PRODUCED WHEAT TORTILLAS

Wheat tortillas are produced by three methods: hot-press, die-cut, and hand-stretch. Different types of flours are required for different tortilla processes. Dough preparation and ingredient utilization vary among different operations. Each operation involves a unique dough-forming procedure that then requires specific flour characteristics.

Consequently, tortillas have different properties and end uses.

Bejosano and Alviola, 2015

Wheat flour is the major and most-relevant ingredient used for manufacturing flour tortillas and **the quality of the finished product depends greatly on the quality of the flour.**

Wheat flour requirements are determined by the desired tortilla characteristics, the formula, processing conditions, and equipment.



# Wheat Flour Tortillas

---

Tortilla flour is considered by milling companies to be a specialty type of flour. It is milled for the tortilla industry as a grade separate from bread, pastry, or all-purpose flours.

Enriched, lightly bleached, hard-wheat flour is generally used for tortillas. It is made from proprietary blends of hard red winter wheat.

Bejosano and Alviola, 2015



Hard Red Winter



Versatile, with excellent milling and baking characteristics for wheat foods like hearth breads, hard rolls, croissants and flat breads.



# Wheat Flour Tortillas

---

## MAJOR FLOUR CONSTITUENTS AFFECTING TORTILLA QUALITY



### Protein Quantity & Quality

Most flour tortillas are made using all-purpose flours containing intermediate protein content of around 11%. The milling operations use wheat blends in order to obtain the most suitable flours for this application.

Flours with stronger protein quality are used in hand-stretch and die-cut tortilla operations, while a wider variety of flours is used in hot-press operations.

A strong-protein flour makes tortillas with smaller diameter but with longer shelf stability.

On the other hand, a weak protein-strength flour makes larger-diameter tortillas that have short shelf stability.

**Thus, flour with intermediate protein quality would be appropriate.**

**Flour protein quality appears to have greater impact on tortilla properties than does protein quantity!**

# Wheat Flour Tortillas

## MAJOR FLOUR CONSTITUENTS AFFECTING TORTILLA QUALITY

### Starch


**Damaged starch**, Starch particle size, Amylose content, Retrogradation

Starch damage in flour affects the properties of tortillas. As starch damage increased, flour tortillas becomes less stretchable, and firmness and rollability increases.

To make a tortilla with acceptable rollability, starch damage should not be too high.

Mao and Flore, 2001

For any product, there is an optimum level of starch damage!



**Starch damage affects water absorption, stickiness, consistency, flexibility and retrogradation!**



# Wheat Flour Tortillas

---

## MAJOR FLOUR CONSTITUENTS AFFECTING TORTILLA QUALITY

### Starch


Damaged starch, **Starch particle size**, **Amylose content**, Retrogradation

Flours with larger starch granules, were found to produce better tortilla texture.

Mao and Flore, 2001

Flour amylose content greatly affected tortilla textural characteristics.

Guo et al., 2003



**Starch affects rollability, appearance (cracks, texture), stickiness, shelf life!**

# Wheat Flour Tortillas


## MAJOR FLOUR CONSTITUENTS AFFECTING TORTILLA QUALITY

### Starch

Damaged starch, Starch particle size, Amylose content, **Retrogradation**

After baking, the starch will tend to partially recrystallize. This phenomenon is called retrogradation and explains why the products become hard (stale).

The faster the starch retrogradation, the faster the tortillas will lose its freshness. As a result, flours with slow retrogradation are favored.



Wang et al., 2015



# Wheat Flour Tortillas

---

## PHYSICAL DOUGH PROPERTIES/DOUGH RHEOLOGY

**The main techniques used for measuring cereal dough properties have traditionally been divided into descriptive empirical techniques and fundamental measurements.**

Within the baking industry there is a long tradition of using descriptive empirical measurements of rheological properties, with instruments such as the Penetrometer, Texturometer, Consistometer, Amylograph, Farinograph, Mixograph, Extensograph, Alveograph, Mixolab various flow viscometers and fermentation recording devices.

They have provided a great deal of information on the quality and performance of doughs such as consistency, hardness, texture, etc.

Dobraszczyk, 2020

# Wheat Flour Tortillas

---

## PHYSICAL DOUGH PROPERTIES/DOUGH RHEOLOGY

**Water Absorption:** Tortillas require less water compared to bread doughs. Water absorption affects; Size & Shape, Appearance (Pillowing), Stickiness (Overhydration), Shelf Life!

**Mixing Time:** Tortillas require less water compared to bread doughs. The under-and overmixing modify dough texture, gluten development, and tortilla quality. Mixing time affects; Size & Shape, Appearance, Stickiness, Shelf Life!

**Extensibility:** Extensibility is the capacity of the dough to be stretched without breaking. Extensibility affects; Size & Shape, Appearance!

# Wheat Flour Tortillas

---

## PHYSICAL DOUGH PROPERTIES/DOUGH RHEOLOGY

**Elasticity:** Elasticity is the tendency of the dough to return to its initial position after its shape is distorted, such as by pressing.

Elasticity affects; Size & Shape, Appearance!

**The rheological properties of dough determine its behavior during dividing, rounding and molding, as well as the quality of the finished products!**

# Wheat Flour Tortillas

---



## PHYSICAL DOUGH PROPERTIES/DOUGH RHEOLOGY

### A Mixolab Study

“The Mixolab profiler showed that a good flour for hot-press tortillas had a relatively lower absorption and short dough mix time compared with a bread flour and should have a significantly higher gluten compared with an all-purpose flour.

Compared with bread flour, the tortilla flour had higher retrogradation and viscosity values.

The Mixolab profiler proved to be a good preliminary test to evaluate flours for hot-press tortillas.

This instrument showed that the best performing flours had a relatively lower absorption and short dough mix time compared with bread flour and a high gluten profile within the category of all-purpose flours.”

Posner et al., 2014

# Wheat Flour Tortillas

## TORTILLA QUALITY

**The ultimate measure for determining the baking performance of flours for tortillas is the baking test.**

- Good quality tortillas are symmetrical, uniform and opaque with toasted spots!
- They should also be soft, flexible without cracking when folded, and puffed!
- Good quality wheat flour tortillas usually have large diameters (17-18 cm) and more than two weeks of shelf stability!



de Barros, 2009

Good-quality tortillas have been defined as tortillas that are symmetrical, uniform, opaque, toast-marked, puffed, soft, flexible without cracking, and having a long shelf life.

Brooker, 2015



Diameter, thickness, weight, specific volume, moisture, rollability, flexibility, appearance (crack/break, blisters), pillowing, color, opacity, firmness, ...!





## Product Inspection



### EXAMPLES **Tortilla, Flat Bread, Pizza Crust**



#### PRODUCTS

**Corn/Flour Tortilla**, Flat Bread, Pizza Crust

#### LINE CONFIGURATIONS

Line Width: 900-2300mm  
Throughput: up to 1500pcs/min

#### SYSTEM MODULES

Top Side Inspection (Color & 3D)  
Bottom Side Inspection (Color)  
Automatic Learning  
Real Time Statistics  
Data Collection & Reporting  
Pneumatic Rejection System  
Hygienic Stainless Steel Frame

#### MEASUREMENT CAPABILITIES

2D/Shape (diameter, roundness)  
3D/Height (peak, mean, slope)  
Bake Color (Lab & BCU)  
Topping/Toast Marks Conformity  
Rim Conformity

#### DEFECT DETECTION CAPABILITIES

Misshapen products  
Edge Defects (bites, straight edge)  
Out of spec (e.g. small, too dark ...)  
Spots (light, dark, foreign mat.)  
Topping defects (e.g. voids)  
Holes

**EyePro System products are leading the effort for incorporating vision technology in the baking and snack food industry.**





## Software Features

ON-LINE measures for each product



### SHAPE/2D CONTROL

- LENGTH / WIDTH
- DIAMETER / AREA
- ASPECT / EDGES
- BITES / TAGS
- HOLES BURNT SPOTS
- OVERLAPPING PRODUCTS
- ROUNDNESS

The screenshot shows the EyePro System software interface. The main window is titled "Display" and shows a model of "TORTILLAS". A central image of a tortilla is displayed with a green bounding box. To the right, a panel titled "Analysis results: product #2" lists various parameters and their values:

Category	Parameter	Value
POSITION:	Position (X, Y center)	(mm) (175, 441)
	Position (X, Y center)	(mm) (527, 407.8)
	Bounding box (X, Y, Z)	(mm) (147, 241)
SHAPE (Area):	Area	(mm²) 33334.4
	Major axis	(mm) 209.6
	Minor axis	(mm) 162
	Aspect ratio	1.352
	Roundness	0.965
	Triang. factor	0.967
	Perimeter of hole	0
	Perimeter of hole	0.967
SHAPE (Perim):	Min. radius	(mm) 85.2
	Max. radius	(mm) 100.3
	Max. Area. radius	70.4
	Distance from object center	(mm) (0.865)
	Radius and Area. difference	3.4
	Max. Area. diff.	882
MEASURES:	Min. diameter	(mm) 101.9
	Max. diameter	(mm) 209.6
COLOUR:	Black colour (mean, BCU)	0.01536
	Colour average	28.9%
	Black difference from mean	0.0
SPOT:	Light spot area	(mm²) 0.0
	Dark spot area	(mm²) 17.9
TOASTED MARKS:	Dark spots	0.0%
	Mark's regularity	0.0%, 0.0%, 0.0%, 0.0%
	Mark's avoidance (X, Y of point)	0.0%

At the bottom of the interface, there are controls for "Thresholds", "Colours", and "Toasted marks", along with a "Zoom" button. The status bar at the bottom shows coordinates and file path: "XY = (1238, 413) RGB = (13, 32, 102) M = 49 C:\CurrentWork\Q6aleTortillas\images\RGB1\tortillas02.bmp".



## Software Features

ON-LINE measures for each product



## FULL COLOR ANALYSIS

- CIE-LAB OR BCU
- BAKE COLOR INTENSITY (MEAN/EDGE/CENTRE)
- TOPPING (COLOR/DISTRIBUTION)
- BURNS / RESIDUES
- SPOTS / HOLES
- SURFACE DEFECTS (TEXTURE)

The screenshot shows the EyePro Learning software interface for tortilla analysis. The main window displays two circular images of tortillas with green and red overlays. The right-hand panel shows analysis results for 'Products inner part (base and toasted marks)' and 'Products outer part'.

Colour	L	a	b	%	Marks
01.50	-2.05	7.56	43.33		
00.03	-3.98	11.67	23.12		
70.17	-3.14	8.40	14.65		
75.80	-1.74	17.45	6.59		X
75.30	-2.75	6.93	4.25		
75.06	-2.52	12.16	3.69		
64.87	3.07	20.41	3.04		X
69.70	-1.59	9.51	1.30		X

Colour	L	a	b	%	Marks
80.25	-3.17	8.93	89.04		(base)
72.08	-0.38	17.33	10.83		(toasted marks)
-	-	-	0.04		(spots)

**Products outer part**

**Edge**


Edge width (mm):

Colour	L	a	b	%	Marks
70.56	-2.97	8.59	-		(edge)

# Software Features

## Defect detection and automatic rejection





**Analysis results: product #1**

**POSITION:**  
 Position [X,Y centre] (pixel) [274, 410]  
 Position [X,Y centre] (inches) [9.1, 13.5]  
 Bounding box [dX,dY] (pixel) [281, 299]

**SHAPE:**  
 Area (inches<sup>2</sup>) 71.7  
 Major axis (inches) 9.9  
 Minor axis (inches) 9.3  
 Aspect ratio 1.067  
 Eccentricity 0.937  
 Shape factor 0.930  
 Number of holes 1  
 Shape class shape  
 Likeness 99%

**EDGE:**  
 Min. radius (inches) 4.6  
 Max. radius (inches) 5.1  
 Min./Max. radius 89.1  
 Distance from object centre (inches) 0.031  
 Radius vect max. difference 0.4  
 Likeness 99%  
 Straight edge presence absent  
 Maximum irregularity area (inches<sup>2</sup>) 0.0000

**MEASURES:**  
 Min. diameter (inches) 9.1  
 Max. diameter (inches) 9.8  
 Mean diameter (inches) 9.5  
 Roundness 93  
 Hole perimeter 2.0

**COLOUR:**  
 Bake colour (mean, BCU) 72.9 4.85  
 Colour evenness 97.5 %  
 Bake difference from mean 0.0


**SPOT:**  
 Light spot area (inches<sup>2</sup>) 0.000  
 Dark spot area (inches<sup>2</sup>) 0.000  
 Foreign material area (inches<sup>2</sup>) 0.000

**TOASTED MARKS:**  
 Marks quantity 4.8 %  
 Marks repartition 5.5 %, 4.6 %, 5.3 %, 3.6 %  
 Marks avoidance (% of zones) 0.0 %

← HOLE DEFECT

# Software Features

## Defect detection and automatic rejection

**Analysis results: product #2**

**POSITION:**  
 Position [X,Y centre] (pixel) (1769, 416)  
 Position [X,Y centre] (inches) (58.5, 13.7)  
 Bounding box [dX,dY] (pixel) (239, 271)

**SHAPE:**  
 Area (inches<sup>2</sup>) 54.8  
 Major axis (inches) 10.8  
 Minor axis (inches) 6.5  
 Aspect ratio 1.658  
 Eccentricity 0.603  
 Shape factor 0.536  
 Number of holes 0  
 Shape class **NO CLASS**  
 Likeness 0%

**EDGE:**  
 Min. radius (inches) 3.1  
 Max. radius (inches) 5.5  
 Min./Max. radius 56.0  
 Distance from object centre (inches) 0.015  
 Radius vect max. difference 0.0  
 Likeness 0%  
 Straight edge presence present (~ 212°)  
 Maximum irregularity area (inches<sup>2</sup>) 0.6830 (315° -> 346°)

**MEASURES:**  
 Min. diameter (inches) 6.1  
 Max. diameter (inches) 10.4  
 Mean diameter (inches) 8.1  
 Roundness 59  
 Hole perimeter 0.0

**COLOUR:**  
 Bake colour (mean, BCU) 0.0 0.00  
 Colour evenness 100.0 %  
 Bake difference from mean 74.0

**SPOT:**  
 Light spot area (inches<sup>2</sup>) 0.000  
 Dark spot area (inches<sup>2</sup>) 0.000  
 Foreign material area (inches<sup>2</sup>) 0.000


**TOASTED MARKS:**  
 Marks quantity 0.0 %  
 Marks repartition 0.0 %, 0.0 %, 0.0 %, 0.0 %  
 Marks avoidance [% of zones] 0.0 %

← **MISSHAPEN PRODUCT**



Vision Technology used to automatically adjust process equipment

Example on Tortilla applications



### Automatic Process Control

Flour tortilla

- Vision System measures product diameters and automatically adjust the press settings
- Improved product size consistency
- Waste reduction

# Vision Inspection Systems

## Tortilla Measurements

Sightline Process Control products are specialized, 3D color vision technologies used in inspection systems and automation solutions for process control.

### Color Analysis

- Avg Color of Top
- Avg Color of Bottom (toast/dark marks ignored)

### Anomaly Detection

- Holes, Folds, Tears, Bites, Tails
- Irregular Edges, Foreign Objects
- Dark Spots & Burn Marks

### Contour Defects




- Edge Roundness
- Blowout Detection
- Bites, Tails, Folds
- Curling

### Toast Marks

- Color Analysis
- Total Number
- Area/Coverage

### 2D Geometry

- Min, Max, Avg Diameter
- Roundness/Shape Verification
- Area Measurements



a KPM ANALYTICS brand

# Corn Tortillas

---

## TYPES OF COMMERCIALLY PRODUCED WHEAT TORTILLAS

### Lime-Cooked Maize Products

Three basic types of products are industrially produced from lime-cooked maize: table or soft tortillas, corn chips, and tortilla chips.

Corn and tortilla chips are primarily produced and consumed in developed countries, where they have an important share of the salted snack-food market.

Modern production of tortilla chips and corn chips has evolved into specialized processes. They are quite different from the original tostadas made by frying stale, maize table tortillas.

Rooney and Serna-Saldivar, 2016



# Corn Tortillas

---

## MAJOR CORN CONSTITUENTS AFFECTING TORTILLA QUALITY

### The types of corn

Corn is classed according to color and kernel configuration into yellow or white and dent or flinty.

White, yellow, and blue corns are processed into tortillas and chips in the United States and Mexico.

The main criteria used to select corn for the manufacturing of tortillas and related snacks are related to their physical properties because they affect chemical composition, functionality, processing parameters, and end use. The preferred corns should be medium to large sized (290–340 g/1,000 kernels) and possess a round crown, a smooth dent, an intermediate endosperm texture, an easily removable pericarp, a clean bright white or yellow color, and a kernel with tolerance to damage during handling.

**Whole, sound, mature kernels of maize with a high proportion of hard to soft endosperm yield more masa and tortillas after nixtamalization!**



# Corn Tortillas

---

## MAJOR CORN CONSTITUENTS AFFECTING TORTILLA QUALITY

### **The types of corn**

Table tortillas are preferably manufactured from dent white maize, although yellow kernels are also used or blended with white kernels before cooking. The grain's physical properties play the main role in the quality of table tortillas and related products because they affect processing parameters and influence end-product features.

The main quality control parameters for fresh masa and dry flours are particle size distribution, water absorption, water solubility, viscosity, color, and pH. Likewise, the principal quality measurements for wheat flour tortillas are dough rheological properties determined with the Farinograph, Alveograph, Mixolab, and related instruments.


Brooker, 2015

# Corn Tortillas

## MAJOR CORN CONSTITUENTS AFFECTING TORTILLA QUALITY

### Assessment of the quality of fresh nixtamalized maize doughs

Mixolab discerns among different degrees of cooking in nixtamalized maize doughs (NMD),  
The milling degree did not affect the Mixolab profile,  
Mixolab parameters were strongly correlated with starch damage in NMD,  
Mixolab showed more correlations with starch and tortilla features compared to RVA,  
Mixolab resulted an adequate method to evaluate the quality of fresh NMD!



Johanan et al., 2021

# Corn Tortillas

## CASE STUDY – CORN TORTILLA

Partnership with a tortilla manufacturer:

- The partner is making tortilla from nixtamalized corn masa,
- Difficulties are often faced during the process,
- They appear to be linked to the corn masa's quality,
- The Mixolab can be used with **sampled doughs** thanks to the dough kit

3 samples of various quality (different milling procedures) were tested

Objectives :

Check if the Mixolab is able to analyze corn masa.  
Is there a signal?, Is it repeatable?, Is it discriminant?

If yes, is the Mixolab able to classify samples according to their qualities?

Bosc-Passive, 2021

1 - Good



2- Medium




3- Bad



# Corn Tortillas

## CASE STUDY – CORN TORTILLA

- ✓ A clear signal is detected, the curves are complete
- ✓ The repeatability curves are well overlaid
- ✓ In general, the lower the torque, the better the masa.



Torques Cs, C3 and C4 allow to classify efficiently the different qualities of doughs in a discriminative way.

**When the dough is heated, the lower the torques, the better the sample is according to the process performances**

**Discrimination OK !  
Feasibility confirmed!**

Bosc-Passive, 2021

# Tortilla Quality

	NIR	SD MATIC	SRC-CHOPIN	ALVEOLAB	MIXOLAB
<b>Size &amp; shape</b>			<b>++</b> absorption	<b>+++++</b> Extensibility/Elasticity	<b>+++</b> Hydration
<b>Rollability</b>	<b>+</b> content	<b>+++</b> High level	<b>++</b>	<b>+++++</b> Protein quality	<b>++++</b> Amylase, starch damage
<b>Appearance</b>			<b>+++</b> Damaged starch	<b>+++++</b> Protein on crack, edges, pillowing, exture	<b>+++++</b> Starch properties action on Cracks and texture
<b>Stickiness</b>	<b>+</b> content	<b>++</b> High level	<b>++</b> High level	<b>+++</b> Protein quality	<b>++++</b> Starch properties, hydration
<b>Shelf Life</b>	<b>+</b> content	<b>+++</b> High		<b>++++</b> Protein quality	<b>+++++</b> Starch properties and damaged starch
<b>Ingredient evaluation</b>				<b>+++++</b> Ready	<b>+++++</b> Ready
<b>Work on complete formulas</b>				<b>+++++</b> Ready	<b>+++++</b> Ready
<b>Work on dough coming from production line</b>				<b>++</b> Possible	<b>+++++</b> Ready

# Conclusion

---



While the tortilla quality depends on other processing variables and the formulation that the manufacturer uses, since wheat flour is the significant and most relevant ingredient used for manufacturing flour tortillas, its quality dramatically affects the final product quality.

Likewise, for corn tortillas, the quality of maize for alkaline cooking is critically important to produce high-quality products.



# Conclusion

---

The objective methods are reliable, sensitive, and time-saving and replace subjective measurements.

In the baking industry, various objective methods are used to characterize the rheology of wheat flour dough and final products. These objective measurements characterize the rheology of wheat flour dough and tortillas.

**There is a need to develop standard approved predictive quality tests for tortillas! In addition, there is a need to develop wheats with improved functionality for the tortilla and flatbread market!**





# Global Leader in Analytical Instrumentation

*KPM Analytics enables companies in agriculture, food, and related industries to effectively manage quality and protect value of their brands*

## Summary Highlights

- Established in 2015, HQ in Boston, MA, USA
- 7 strong brands located in four countries
- Long term partnership with blue chip customers in various industries around the world
- Wide product range for R&D, quality control and process optimization
- ~50 patents granted or pending, ~40 trademarks
- Global applications experts with decades of experiences
- Over 15,000 Installs in 100+ countries
- 14 direct sales & service offices globally
- Over 200 sales and service distribution partners worldwide

## Industries



Agriculture



Food & Beverage



Feed, Forage &  
Animal food



Environmental  
Testing



Clinical  
Diagnostics



Industrial  
Manufacturing


## KPM Brands







# Timeline and Progression of KPM Analytics

KPM Analytics brands have long and successful histories servicing customers worldwide



# Impacting Quality Across the Entire Value Chain



From raw ingredients to prepared foods, our products rapidly assess critical properties and overall quality




# Trusted Partner Supporting Every Step of the Food Production Process

KPM brand solutions are used in every process step throughout the grain, milling and baking supply chain


Breeding  
Through Elevator




Milling




Baking




# CHOPIN Technologies is a brand of KPM Analytics



INSPECTION ▶ 2D-Shape ▶ Color ▶ 3D-Height ▶ Top & Bottom ▶ REJECTION



# Extensive Global Sales and Service Network



THANK YOU FOR YOUR ATTENTION  
ANY QUESTIONS ?

---

[hboyacioglu@kpmanalytics.com](mailto:hboyacioglu@kpmanalytics.com)

