



Tortilla Wheat Flour characteristics and quality

A. Dubat, Business development Director

adubat@chopin.fr

September 15, 2018





Who are we?



NIR process moisture measurement

- Animal food
- Paper converting
- Chemicals & Pharmaceuticals
- Human food
- Tobacco
- Wood
- Minerals
- Textiles
- Renewable energy

NIR offline analyzers

- Agricultural
 - Feed ingredients
 - Compound feed
 - Pet food
 - Oilseeds
 - Forage
- Food and Dairy
- Environmental & Industrial
 - Water
 - Soils

Methods and equipment for the analysis of characteristics of cereals and derivatives

- Sample preparation
- Compositional analysis
- Functional analysis
- Bin temperature monitoring

Instrumentation, reagents and tests

- Medical diagnostics
- Environment
- Food & beverage industries





Our Customers

CHOPIN Technologies customers come from across the cereal chain, from seed to finished product

- Breeders
- Traders
- Elevators
- Primary processing industries
 - Milling and malting
- Manufacturers of ingredients or additives for flour
- Baking industries, secondary processors
 - Bakers
 - Pasta producers
 - Snack producers
- Laboratories, Research Institutes, Universities

CHOPIN Technologies DNA: Innovation



Mixolab

Measures the characteristics of dough during mixing, heating and cooling determining the quality of starch and protein



Quatuor II

Automatic dockage tester for determination of seed impurities



Aquaneo

Moisture meters for the certified measurement of the moisture content of cereals and oilseeds



Mixolab 2

New generation of Mixolab

Alveolab

New generation of Alveograph, full version

AlveoPC

New generation of Alveograph, basic version

SRC CHOPIN

Automated measurement of the solvent retention capacity

2003

2005

2007

2009

2010

2011

2013

2014

2016

2017

2018

SDmatic

Analyzer of damaged starch in less than 10 minutes



Infraneo

Infrared analyzer for whole grains and powdery products



Infraneo JR

Infrared analyzer with a smaller capacity for whole grains and powdery products



RheoF4

Unique solution to analyze dough proofing properties



LabMill

Reproduces industrial mill productivity. Analyzes soft and hard wheat milling capacity



AmyLab FN

Faster method for prediction of Hagberg falling number value



Spectralab

High performance NIR analyzer designed for process control in flour mills



Quality: Our Priority

- CHOPIN Technologies operates a Quality Management System which complies with the requirements of ISO 9001 : 2015 standard.
- CHOPIN Technologies manufactures high quality equipment, recognized by numerous international organizations.
- CHOPIN Technologies takes part, as an expert, in different standardization groups.




GOST



GB/T



Global presence, both direct and through distributors



September 15, 2018

Let's talk about Tortillas and Quality



Method

Basis of the work supporting this presentation

- 1/ Our experience on flour quality and how to measure it
- 2/ Analysis of previous TIA Technical Conferences presentations
 - Steve Bright: 2017 (2), 2015 (2) & 2016
 - Cristina Primo Martin: (2017)
 - Tom Jondiko: (2016)
 - Dilek Austin: (2016)





Main flour quality - related identified control point

- Tortilla Size and Shape
- Rollability/Foldability/Flexibility
- Appearance for users
 - Cracks
 - Edges
 - Pillowing
 - Texture
- Stickiness
- Shelf Life

Tortilla Size and Shape

Tortilla too small

- Strong flour : **Elastic**
- Under mixed : **Elastic**
- Underhydrated : **Elastic**
- (Cold dough, under scaling, excessive floor time after mixing, poor press setup, oven shrinkage)


Tortilla too large

- Overly **extensible** flour
- Flour quality
- Protein quantity/quality
- **Overhydration**
- (Over mixing, hot dough, high level of reducing agents, press too severe, high fat level, over scaling)



Tortilla Size and Shape

- Importance of Gluten Quality
 - Strong flour tend to produce tortillas with smaller diameter
 - Weak flour can give;
 - Thinner tortilla, less **gas retention**
 - Severe Sticking
 - Fragile Tortilla, poor shelf stability
 - Poor shape (too large)



Sources : Steve Bright TIA technical conference 2015
Cristina Primo Martin TIA technical conference 2017


Tortilla Size and Shape

- By modifying **gluten properties**, proteases or oxydases are enzymes that help control Tortilla diameter.




40 ppm
L-Cysteine

Tortilla dough
conditioner




Rollability/Foldability/Flexibility

- Protein is important...
- Low Protein Content (<9%)
 - Tortilla cracks easily
 - Larger diameters
- Protein content >12%
 - Better foldability
 - But too long mixing & resting time
 - Smaller diameters




Rollability/Foldability/Flexibility

- ...and Starch is also playing a major role



Rollability/Foldability/Flexibility

- ...as well as starch quality and damage



September 15, 2018

Flour: Starch Quality

- Good quality tortillas produced using wheat flour of intermediate protein and low level of starch damage (Waniska et al 2004) .

Starch characteristics:

- Measured using rapid visco analyzer.
- Peak and final viscosity.
- Affect cooking conditions.
- Tortilla flexibility.
- Tortilla shelf stability.




SOLVAIRA SPECIALTIES

CHOPIN
TECHNOLOGIES

Appearance for users

CRACKS

- ◆ Low protein content (< 9%):
 - gives tortillas that crack easily,



*Enzyme : amylase

EDGES

Reducing agents

Used to aid pressing the dough to a desired size

- L-Cysteine
- Sodium metabisulfite
- Inactive yeast

Low use rates (ppm) cost effective.

Under use

- Small sizes
- Brittle, rough edges
- Laced edges
- \$\$\$ Cost in rejects

Sources : Cristina Primo Martin TIA technical conference 2017
Dilek Austin TIA Technical Conference 2016
Steve Bright TIA Technical Conference 2017

Appearance for users

TEXTURE

PILLOWING

Strong Flour
enhances pillowing
gas retention



- ◆ Texture properties can be affected by modifying starch structure: A limited starch breakdown is needed to give the tortillas the required strength and eating properties (bite, chewiness)

Mouthfeel, bite

- ✓ Short tender bite
 - Established by formula and process
 - Lamination*
 - From leavening
 - Not over pressed –pressure, dwell time, temperature
- X Leathery, tough bite
 - High translucency*
 - Insufficient leavening
 - Hot press
 - Extended press dwell times



Sources : Steve Bright TIA Technical conference 2015
Cristina Primo Martin TIA technical conference 2017

Stickiness



- Dough Stickiness
 - Starch Damage, Protein quality

◆ Damaged starch: starch granules can be damaged during milling. They have a very high water absorption that can lead to stickiness and quality problems.

- Tortilla stickiness
 - Protein quality
 - Formulation (fat, sugar...)

Sticking - Reducing Agents

L-Cysteine and sodium metabisulfite

- greater extensibility in the dough
- higher levels (>60ppm) lead to weak protein and crust resilience.
- Increases the occurrence of sticking

- Obtain dough consistency through
 - full mix development
 - Optimizing flour to water ratio

Shelf Life

Wheat Flour - Starch

Changes occurring in the starch during baking affect their functionality:

- ◆ After gelatinization starch tends to regain crystalline structure: Retrogradation
- ◆ Negative effect on shelf life (staling): Fresh tortillas have good foldability but during storage, tortillas get harder and loss foldability.
 - ◆ Amylose retrogradates fast after baking (\approx 30-60 min)
 - ◆ Amylopectin retrogradation starts one day after baking

Preservative / Shelf Life Extenders

PRESERVATIVES - Low use rate

- Short shelf life
- Moldy tortillas

Gums

Added for shelf life extension and anti-sticking

- Many types of gums available
- Guar, Cellulose, Xanthan

Emulsifiers

Added for shelf life extension and anti-sticking

- Many types of emulsifiers available
- Monoglycerides, SSL, DATEM, Lecithin

Insufficient gluten (in flour and added as an ingredient)

- Thinner tortillas, less gas retention
- Severe sticking
- Fragile tortillas, poor shelf stability. No elasticity in the finished tortilla
- Poor shapes in the press, excessive sizes


Sources : Steve Bright TIA Technical conference 2015
Cristina Primo Martin TIA technical conference 2017

Summary

Flour-related main control points

	Protein	Starch	Water
Size & shape	++++ (extensibility/Elasticity)		++ (Overhydration)
Rollability	+++ (Protein content...and quality)	++++ (Properties improved with Amylase action) +++ (Starch damage)	
Appearance	++++ (Cracks, Edges, Pillowing, Texture)	+++ (Cracks, texture)	++ (Pillowing)
Stickiness	++++ (Protein content...and quality)	+++ (Starch properties)	++ (Overhydration)
Shelf Life	++++ (Protein content...and quality)	+++ (Starch Properties) +++ (Starch damage)	

September 15, 2018



Main flour component impacting Flour quality for Tortilla Production

- **PROTEIN**
 - Content, Strength, Elasticity, Extensibility
- **STARCH**
 - Properties (viscosity, retrogradation), starch damage
- **HYDRATION**
 - Protein content, starch damage
- **FORMULA**
 - How to test complete formula, can we test dough from the processing line?

Best potential tools



Infraneo



SDmatic



Alveograph



Mixolab



Tortilla flour quality

Quantitative Analysis: NIRS





Reminder on NIR basics (1/2)

- NIR has many advantages :



- Simplicity
- Speed
- Clean
- Performances...

- But keep in mind these devices are measuring some light wavelength absorbed that depends on the vibration properties of some known chemical linkages (ex : NH for the protein).



Reminder on NIR basics (2/2)

- The absorption peaks are mathematically analysed according to the sample concentration for the requested element to analyse (ex : proteins)
- Important consequences:
 - The infrared method relies on the accuracy of the reference method which is mirrored
 - The more a parameter or a component is far for the regular NIR conditions, the more its estimation becomes difficult. The calibration developed for such components are giving « trends » results and can't be that accurate.



NIRS for Flour

	Specific absorption peaks ?	Relations with a component with a dedicated peak ?	Comments
Moisture	Yes (many)	No	Measurement
Proteins	Yes (many)	No	Measurement
Wet Gluten	No	Direct (proteins)	Measurement
Water absorption capacity	No	Indirect (Proteins, starch)	Estimation
Ashes	No	indirecte (relation with Cellulose)	Estimation
Starch	Yes	No	Measurement


BEWARE of « calibrations » not based on sound science, these can be strongly misleading (Alveograph, bread volume... data for example)

Measuring Starch Damage



1. Introduction - starch

- Starch is semi crystalline structure composed of two polymers : amylose and amylopectin
 - The grain kernel contains 65-70% starch
 - Starch is the main element of flour; from 68 to 72%




2. Starch Damage

- During milling, the starch granules are more or less damaged depending on:
 - Wheat hardness (**genetic criteria**)
 - Milling adjustments , including wheat tempering (**mechanical criteria**)



3. Why measure starch damage?

- For any product, there is an optimum level of starch damage.



3. Why measure starch damage?

- A lack of control of starch damage may lead to several issues during product manufacturing.



Stickiness



Fermentation



Texture



Color



Cracks

3.1 Effects on water absorption

- Protein absorbs 1.8 times its weight of water.
- Pentosans :10 times
- Native starch : 0.4 times
- Damaged starch: 3 to 4 times



3.2 Effects on biscuits cracks


- An excess of starch damage causes:
 - Broken biscuits when packaging is opened!
 - Biscuit with colour defaults (too dark / too light)
 - Non standard size





Dough Rheology : The Alveolab





The Alveograph is all about Bubble!



- 1- system in position
- 2- blowing starts
- 3- expansion
- 4- until bubble bursts



Flour Strength (W)




Baking strength depends on :

- Protein quality and quantity
- Starch damage
- Enzymatic activity
- And more

-It is a « global view » it is always better to work on P, G and le.

Strength is composed of different aspects




Elasticity (Ie)



- Elasticity is the hability of a dough to return to its original state when the stress disappears.
- It is measured on the alveograph 40 mm after the test beginning.
- At this time, trhe instrument has injected 200 ml of air inside the bubble.
- The volume of the bubble at this very moment indicates dough resistance to deformation which is a measure of elasticity

Focus on elasticity measurement




Figures given as example ; do not correspond to real values!

ICC Pretoria February 2010



Processing aids action : example on Cysteine




Dough Rheology: The Mixolab



What is a Mixolab


The Mixolab measures the characteristics of dough during mixing (water absorption, development time, stability, etc..) while evaluating the quality of protein and starch.



User benefits:

By heating up and cooling down the dough in one test, you have a very complete dough rheology analysis.


Standard Mixolab curve




Compliant with

ICC 173/1 ; AACC 54-60-01 ; NF V03-765 ; ISO 17718:2013 ; GOST R 54498-2011

Mixolab Profiler : make it simple



CHOPIN




Target
Min. index
5-55-543
Max. index
7-76-664
Mixolab
5-68-686

Index type	Values	Significance : the higher the index value the...
ABSORPTION	From 0 to 9	... more the flour absorbs water
MIXING		... more the flour is stable at kneading
GLUTEN+		... more the gluten resists heat
VISCOSITY		... greater the dough's viscosity when heated
AMYLASE		... weaker the amylase activity
RETROGRADATION		... shorter the cooked product's shelf life will be

CHOPIN
TECHNOLOGIES

Example 1 : α -amylase testing


Fungal



Only flour


Flour + 1 dose

Bacterial




Flour + 2 doses

Maltogenic




Flour + 4 doses


Protein




Gelatinization




Liquefaction



Rétrogradation



Example 2: Gums testing



0,0 % carrageenan


0,5 % carrageenan

1,0 % carrageenan

"Thickening"
Carrageenan



"Gelling"
Carrageenan



And why not testing a dough piece from the production line?


- In this example, dough was taken from the industrial mixer while they were looking at hydration impact on dough texture

- **Industrial mixer:**

- Flour : 80 Kg
 - Water : 36,2 L
 - Water : 36,4 L
 - Water : 36,8 L

- **Mixolab:**

- Dough sample: 100 g
- Rotation speed: 100 rpm
- Mixer temperature: 26°C





Reminder :

Flour-related main control points

	Protein	Starch	Water
Size & shape	++++ (extensibility/Elasticity)		++ (Overhydration)
Rollability	+++ (Protein content...and quality)	++++ (Properties improved with Amylase action) +++ (Starch damage)	
Appearance	++++ (Cracks, Edges, Pillowing, Texture)	+++ (Cracks, texture)	++ (Pillowing)
Stickiness	++++ (Protein content...and quality)	+++ (Starch properties)	++ (Overhydration)
Shelf Life	++++ (Protein content...and quality)	+++ (Starch Properties) +++ (Starch damage)	

September 15, 2018

Take home message:

Which device for controlling flour quality for Tortilla production?

Wheat Tortilla parameter	NIRS	SDMATIC	Alveolab	Mixolab
Size & shape			+++++ Extensibility/Elasticity	+++ Hydration
Rollability	+ content	+++ High level	+++++ Protein quality	++++ Amylase, starch damage
Appearance			+++++ Protein on crack, edges, pillowing, exture	+++++ Starch properties action on Cracks and texture
Stickiness	+ content	++ High level	+++ Protein quality	++++ Starch properties, hydration
Shelf Life	+ content	+++ High	+++++ Protein quality	+++++ Starch properties and damaged starch
Ingredient evaluation			+++++ Ready	+++++ Ready
Work on complete formulas			+++++ Ready	+++++ Ready
Work on dough coming from production line			++ Possible	+++++ Ready

THANK YOU

GRACIAS
ARIGATO
SHUKURIA
JUSPAXAR
DANKSCHEEN
TASHAKKUR ATU
YAQHANYELAY
SUKSAMA
EKHMET
BIYAN
SHUKRIA
TINGKI
GRAZIE
MEHRBANI
PALDIES
BOLZIN
MERCY
GOZAIMASHITA
EFCHARISTO
FAKAAUE
KOMAPSUMNIDA
MAAKE
LAH
MERASTAWHY
GAEJTHO
AGUYJE
FAKAAUE
TAVTAPUCH
MEDAWAGSE
BAUKA
JUSPAXAR
SPASSIBO
SNACHALHUYA
HURUN
CHALTU
WADEEJA
MAITEKA
YUSPAGARATAM
HUI
UNALCHEESH
NENACHALHYA
DENKAUJA
MERSE
SPASIBO
NATUR
GUT
EKOJU
SIKOMO
MAKETAI
MINMONCHAR



www.chopin.fr

CHOPIN
TECHNOLOGIES