



 deutscheback
Together. Better Baking.

Effect of Flour Storage Temperature on Tortillas



Nicolas Charalampidis

R&D



Rheology Farinograph

 STORK TECHNOLOGY CENTER
Futuremiker

Farinograph

- Mixing behaviour
- Water absorption
- Dough development time
- Stability
- Softening



Principle Farinograph





Optimum
consistency

300 g


Interpreting the curve

Farinograph



Interpreting the curve


Farinograph



	Low	Normal	High
Water absorption (%)	<54	56	>60
Dough Development Time (min)	<1.5	1.5-5	>5
Stability (min)	<3	5-7	>10
Softening (BU)	<30	30-70	>70


Farinograph

22°C




Farinograph

2°C storage




Farinograph -18°C storage



Farinograph

-18°C + 3hrs room




Results

	Water absorption Corrected 500FU	Development time	Stability
Flour room 22°C	53.5 %	2.3 min	7.6 min
Flour 2°C	55.5 %	2.8 min	11.7 min
Flour -18 °C	57.8 %	2.4 min	10.2 min
Flour -18 + 3 hours room	54.8 %	2.8 min	10.3 min


Farinograph

22°C




Farinograph

2°C storage




Farinograph -18°C storage



Farinograph

-18°C + 3hrs room



Results


	Water absorption Corrected 500FU	Development time	Stability
Flour room 22°C	57.5 %	1.9 min	4.5 min
Flour 2°C	60 %	2.5 min	11 min
Flour -18 °C	60.7 %	2.2 min	9.3 min
Flour -18 + 3 hours room	59.8 %	2.7 min	13.6 min



Rheology Farinograph +


Farinograph

22°C




Farinograph

2°C storage




Farinograph -18°C storage



Farinograph

-18°C + 3hrs room




Results

	Water absorption Corrected 500FU	Development time	Stability
Flour room 22°C	59.2%	2.4 min	9.6 min
Flour 2°C	60.9 %	2.7 min	10.9 min
Flour -18 °C	61.2 %	2.5 min	8.3 min
Flour -18 + 3 hours room	61.0 %	2.5 min	10.8 min


Farinograph

22°C




Farinograph

2°C storage




Farinograph -18°C storage



Farinograph

-18°C + 3hrs room




Results

	Water absorption Corrected 500FU	Development time	Stability
Flour room 22°C	58.6 %	2.3 min	8.1 min
Flour 2°C	61.2 %	2.2 min	12.7 min
Flour -18 °C	60.9 %	2.8 min	10.1 min
Flour -18 + 3 hours room	61.3 %	3.3 min	13.6 min


Farinograph

22°C




Farinograph

2°C storage




Farinograph -18°C storage



Farinograph

-18°C + 3hrs room




Results

	Water absorption Corrected 500FU	Development time	Stability
Flour room 22°C	59.6 %	2.7 min	9.8 min
Flour 2°C	60.8 %	2.5 min	9.9 min
Flour -18 °C	59.9 %	2.9 min	8.1 min
Flour -18 + 3 hours room	61.2 %	2.7 min	8.5 min


Farinograph

22°C




Farinograph

2°C storage




Farinograph -18°C storage



Farinograph

-18°C + 3hrs room



Results


	Water absorption Corrected 500FU	Development time	Stability
Flour room 22°C	59.9 %	2.4 min	12.0 min
Flour 2°C	60.6 %	2.5 min	9.6 min
Flour -18 °C	60.6 %	2.0 min	8.9 min
Flour -18 + 3 hours room	60.8 %	2.5 min	10.5 min



Rheology Falling number


 STECK TECHNOLOGY CENTER
Futuremiker

Principle Falling number



Interpreting the values

Falling Number



Results

	Falling number	Falling number
Flour RoomTemp	541	465
Flour 2°C	459	446
Flour -18 °C	487	410
Flour -18 + 3 hours room	280	268



Rheology Viscograph


 STEIN TECHNOLOGY CENTER
Futuremiker

Micro-Visco Amylograph

- Flour suspension
- Heating and cooling
- Retrogradation – freshkeeping



Principle MVAG




Measured
Temperature
Resistance
Time

Curve


Interpreting the curve

MVAG



Micro-Visco Amylograph

Flour normal




Evaluation

Point	Name	Time [HH:MM:SS]	Torque [BU]	Temperature [°C]
A	Beginning of gelatinization	00:05:08	10	68.3
B	Maximum viscosity	00:12:18	217	90.0
C	Start of holding period	00:08:00	72	89.5
D	Start of cooling period	00:13:00	216	90.0
E	End of cooling period	00:18:20	241	54.7
F	End of final holding period	00:23:20	254	49.2
B-D	Breakdown		1	
E-D	Setback		25	

Micro-Visco Amylograph

2°C storage




Evaluation

Point	Name	Time [HH:MM:SS]	Torque [BU]	Temperature [°C]
A	Beginning of gelatinization	00:05:18	10	69.4
B	Maximum viscosity	00:11:52	208	90.0
C	Start of holding period	00:08:00	71	89.5
D	Start of cooling period	00:13:00	208	90.0
E	End of cooling period	00:18:20	229	55.0
F	End of final holding period	00:23:20	246	49.1
B-D	Breakdown		0	
E-D	Setback		21	

Micro-Visco Amylograph

-18°C storage




Evaluation

Point	Name	Time [HH:MM:SS]	Torque [BU]	Temperature [°C]
A	Beginning of gelatinization	00:05:40	10	71.9
B	Maximum viscosity	00:11:46	205	90.1
C	Start of holding period	00:08:00	70	89.6
D	Start of cooling period	00:13:00	204	90.0
E	End of cooling period	00:18:20	224	55.3
F	End of final holding period	00:23:20	243	48.8
B-D	Breakdown		1	
E-D	Setback		20	

Micro-Visco Amylograph

-18°C + 3hrs room



Evaluation

Point	Name	Time [HH:MM:SS]	Torque [BU]	Temperature [°C]
A	Beginning of gelatinization	00:05:04	9	67.9
B	Maximum viscosity	00:11:48	207	90.0
C	Start of holding period	00:08:00	72	89.5
D	Start of cooling period	00:13:00	205	90.0
E	End of cooling period	00:18:20	227	55.1
F	End of final holding period	00:23:20	244	49.1
B-D	Breakdown		2	
E-D	Setback		22	



Rheology Gluten

 STORA ENSO TECHNOLOGY CENTER
Futurebakker

Gluten Index

Mixer




Centrifuge



Dryer




Principle Gluten Index



Interpreting the curve

Gluten Index



Sieve

Results

	Wet Gluten	Gluten Index	Dry Gluten
Flour room 22°C	31.1 %	88	10.2 %
Flour 2°C	30.8 %	85	10.4 %
Flour -18 °C	31.4 %	90	10.5 %
Flour -18 + 3 hours room	31.1 %	90	10.4 %



Baking

Tortilla

Ingredients	%
Flour	100
TopBake TFX 100	0.1
Salt	2
Sugar	0.5
Acids	0.4
Calcium Propionate	0.8
Baking powder	2
Shorting	12
Water	58

Production	
Mixing	2 + 12 min
Dough Temp	26°C - 27 °C
Resting time	12 min
Press Temp	210°C
Press time	1.4 Sec
Pressure	700 PSI

Tortilla

+ 22°C

Characteristics	
Dough Temperature	26°C
Dough after mixing	Soft, elastic and not sticky
Dough before pressing	Elastic, soft, a little wet
Comments	Good dough, good press characteristics



Tortilla

+ 2°C

Characteristics	
Dough Temperature	26.2°C
Dough after mixing	Soft, slightly sticky and little bucky
Dough before pressing	Overall ok, a little bucky
Comments	Decent press characteristics, slightly uneven



Tortilla

- 18°C

Characteristics	
Dough Temperature	25.8°C
Dough after mixing	Sticky, tough and bucky
Dough before pressing	Bucky and not relaxed
Comments	Very uneven form and slightly yellow




Tortilla

- 18°C +3H

Characteristics	
Dough Temperature	26.4°C
Dough after mixing	Short, little sticky
Dough before pressing	Bucky
Comments	Uneven press More pillowing



Tortilla





 deutscheback

Together. Better Baking.

Thank you very much for your attention
Nicolas Charalampidis

R&D

ncharalampidis@sterningredients.com.mx
cel. +525559092248